Patient Leaflet

CS - PŘÍBALOVÁ INFORMACE PRO PACIENTA	02
DA - PATIENTINFORMATIONSBROCHURE	05
DE - PATIENTENINFORMATION	08
EL - ΕΝΗΜΕΡΩΤΙΚΟ ΦΥΛΛΑΔΙΟ ΓΙΑ ΤΟΝ ΑΣΘΕΝΗ	11
EN - PATIENT INFORMATION LEAFLET	14
ES - FOLLETO INFORMATIVO PARA EL PACIENTE	17
ET - PATSIENDI INFOLEHT	20
FI - KUVITETTU TIEDOTE POTILAALLE	23
FR - NOTICE D'UTILISATION DESTINÉE AU PATIENT	26
HR - UPUTA O LIJEKU ZA PACIJENTA	29
HU - BETEGTÁJÉKOZTATÓ	32
IT - FOGLIO ILLUSTRATIVO PER IL PAZIENTE	35
NL - BIJSLUITER VOOR DE PATIËNT	38
PL - ULOTKA DLA PACJENTA	41
PT - FOLHETO INFORMATIVO PARA O PACIENTE	44
RO - PROSPECT PENTRU PACIENT	47
SK - PŘÍBALOVÁ INFORMÁCIA PRE PACIENTA	50
SL - PRILOGA Z NAVODILI ZA PACIENTA	53
SR - ИНФОРМАЦИОНИ ЛИСТ ЗА ПАЦИЈЕНТА	56
SV - PATIENTINFORMATION	59

SURGERY

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VLASTNOSTI

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PATIENT INFORMATION LEAFLET

FEATURES

Glubran[®] 2 is a class III surgical medical device (internal and external surgical use) that complies with current European Regulations. Glubran[®] 2 is a liquid of a synthetic nature with a cyanoacrylate base modified by the addition of a monomer synthesised by the manufacturer.

Glubran[®] 2 has pronounced adhesive, sealing and haemostatic properties; once cured, it provides an effective antiseptic barrier against the most common infectious or pathogenic agents in surgery. It is a ready-to-use, transparent, light-yellow liquid. On contact with living tissue and in a humid environment, it cures rapidly, creating a thin elastic film with high tensile strength that ensures solid tissue adhesion. Such a film naturally conforms to the anatomy of the tissues to which it is affixed, is impermeable and is not compromised/lysed/degraded by blood and body fluids. The film can easily be pierced by a suture needle. The curing time varies depending on the type of fabric Glubran[®] 2 comes into contact with, the nature of the liquids present and the amount of product applied. Under conditions of proper application Glubran[®] 2 starts to polymerize after about 1-2 seconds, completing its reaction around 60-90 seconds. At the end of this reaction, Glubran[®] 2 reaches its maximum mechanical strength. Once cured, Glubran[®] 2 no longer has any adhesive power, so that surgical fabrics or gauze can be placed side by side or juxtaposed without the risk of unwanted sticking. In normal surgical procedures, the Glubran[®] 2 film is removed by a process of hydrolytic degradation; the duration of this process varies depending on the type of tissue and the amount of Glubran[®] 2 applied. The polymerisation reaction can generate a temperature of around 45 °C.

Glubran[®] 2 shortens functional recovery time post-surgery and reduces pain treatment of hernias.

Glubran[®] 2 reduces surgery time.

Glubran[®] 2 prevents post-operative complications, such as bleeding, fistula formation of various kinds and collection of biological fluids.

Glubran[®] 2 in endovascular application stops arterial and venous haemorrhages in seconds. Glubran[®] 2 allows the minimally invasive treatment of different types of fistulas.

PACKAGING AND PRODUCT CODES

See the table below. Glubran[®] 2 and its primary packaging are sterile and latex-free.

REF.	Packaging	Format	GMDN	Basic UDI-DI
G-NB2-SUR	10 single-dose/box	1 ml	58777	803190219X
G-NB2-50-SUR	10 single-dose/box	0.5 ml	58777	803190219X
G-NB2S-25-SUR	10 single-dose/box	0.25 ml	58777	803190219X

GLUBRAN[®] 2 IS USED BY QUALIFIED DOCTORS IN PAEDIATRIC AND ADULT PATIENTS:

Glubran[®] 2 must only be used by qualified doctors with experience in using the product.

Glubran[®] 2 has adhesive, sealing, haemostatic, sclerosing, embolising and bacteriostatic effects on tissues. Glubran[®] 2 can be applied to muscle, bone, cartilage, endothelial and adipose tissue, connective tissue, internal epithelia, parenchyma, mucous membranes and serous membranes (peritoneum, pleura, etc.).

It is used in conventional and laparoscopic surgery and in gastrointestinal endoscopy, interventional radiology and vascular neuroradiology procedures. It can be used alone or in combination with sutures, including in patients treated with heparin and patients with hypothermia.

There are no contraindications to the use of this device in paediatric and adult populations.

Please note the warnings and precautions.

Below are some applications of the product. In order to learn the correct way to apply Glubran[®] 2, it is necessary to contact our company or authorised distributor and learn the appropriate technical instructions and professional medical literature.

Cardiac Surgery

- As a haemostatic and suture-reinforcing sealant on the myocardium.
- As a haemostatic and suture-reinforcing sealant on anastomoses (vessel-vascular and prosthetic-vascular).
- As a patch adhesive on the myocardium or to glue the dissection plane in aortic dissections.

Vascular surgery

- As a haemostatic and suture-reinforcing sealant on anastomoses (vessel-vascular and prosthetic-vascular).
- As an adhesive for the treatment of iatrogenic pseudoaneurysms, via percutaneous injections.

Neurosurgery

- As an adhesive and sealant for bonding patches on dural tissue.
- As a sealant to prevent liquor fistulas for trans-sphenoidal closure of the sella turcica.
- As an adhesive and sealant for bonding bone flaps and bone and osteochondral fragments.

Otolaryngology/Maxillofacial Surgery

- As a sealant for the treatment of CSF fistulas in nasal-paranasal and pituitary surgery.
- As a sealant for the treatment of pharyngocutaneous and oroantral fistulas.
- As an adhesive in nose and paranasal sinuses surgery.
- As a sealant of palatal sutures in cleft lip and palate surgery.
- As an adhesive for bonding bone-cartilage fragments.
- As an adhesive to promote the engraftment of skin grafts.

Odontostomatology

- As an adhesive for stabilising the bone window in sinus lift surgery.
- As an adhesive for the stabilisation of bio-filler materials in guided bone regeneration procedures.
- As an adhesive for closing intraoral incisions.

Ophthalmology

- As an adhesive for fixation of the conjunctival flap in pterygioplasty surgery.
- As a sealant and adhesive for closing corneal perforations.

Paediatric surgery

- As a haemostatic on resections or lesions of parenchymatous organs (liver, kidney, pancreas, spleen).
- As a sealant and haemostatic on lung biopsies.
- As a suture substitute in circumcision and hypospadias repair.
- As a sealant and haemostatic in closing organ perforations.
- As a sealant in the treatment of fistulas.
- As a sclerosant in the treatment of gastro-oesophageal varices.
- As an adhesive and sealant in the closure of the patent processus vaginalis.

General surgery

- As an adhesive in the repair of hernias and incisional hernias with patch/mesh in both traditional and laparoscopic surgery.
- As a sealant in the treatment of proctological fistulas.
- As a sealant and reinforcement of both gastrointestinal anastomoses and anastomoses of the terminal rectum.
- As a sealant in the treatment of gastrocutaneous fistulas.
- As a sealant for the prevention of bilomas and biliary fistulas following liver resection.
- As a haemostatic on parenchymal tissues (liver, pancreas, spleen, etc.), on internal epithelia, muscle and adipose tissue.

Bariatric surgery

- As an adhesive for bonding adipose tissue in omentoplasty procedures after sleeve gastrectomy.
- As a sealant and suture reinforcement after Sleeve Gastrectomy or Roux-en Y Gastroenterostomy.

Plastic surgery

- As an adhesive, haemostatic for skin flap adhesion in abdominoplasty procedures.
- As a sealant in the prevention and treatment of seroma in abdominoplasty.
- As an adhesive and haemostatic in rhinoplasty surgery.
- As an adhesive in breast reconstructive surgery.

Thoracic surgery

- As a reinforcing sealant for immediate aerostasis and an improved mechanical seal.
- As a sealant for bronchial and bronchopleural fistulas.

Gynaecological surgery

- As an adhesive and haemostatic in vaginal and perineal plasty.
- As a sealant for the prevention and treatment of inguinal lymphocele and lymphorrhoea.
- As an adhesive to fix meshes for the treatment of prolapse.

Breast surgery

- As a sealant for the prevention and treatment of seroma and lymphorrhoea following mastectomy, quadrantectomy and axillary lymphadenectomy procedures.

Urological surgery

- As a haemostatic and sealant during kidney transplantation.
- As a sealant and haemostatic on lacerations and haemorrhagic renal injuries.
- As a sealant and haemostatic of the excretory tract in partial nephrectomy surgery.
- As a sealant in the treatment of urinary fistulas.
- As a sealant in the treatment of postoperative lymphorrhoea.
- As a suture substitute in phimosis, circumcision and frenulotomy surgery.
- As a graft adhesive in urethroplasty procedures.
- As a haemostatic on endoscopic resections of bladder tumours.

UNDESIRABLE EFFECTS

Although rare, inflammatory reactions may occur at the product application site.

These reactions occur mainly when the amount of Glubran® 2 applied exceeds the recommended dose. In areas prone to infection (such as the urogenital system), an excessive amount of product increases the risk of inflammatory reactions that may cause persistent infections.

In hypersensitive and/or allergic patients, use of the product may cause a severe allergic reaction and, in rare cases, anaphylaxis. In some extravascular applications or during endovascular procedures, sensitive patients may experience a slight burning sensation accompanied by temporary pain, which resolves spontaneously.

If you experience any other adverse effects not mentioned, contact your healthcare professionals.

CONTRAINDICATIONS

Do not apply the product directly to brain tissue or the spinal cord.

Do not apply the product inside a blood vessel, except for treatment during gastrointestinal endoscopy, interventional radiology and vascular neuroradiology procedures.

When used for skin closure, never apply the product within the wound margins, but only to the skin surface, once the thoroughly cleaned wound margins have been perfectly aligned.

Do not use this product in particularly sensitive individuals or pregnant women.

Do not use this product on bleeding varices caused by juvenile liver cirrhosis of unknown origin.

Do not use this product for peripheral nerve anastomoses.

PRECAUTIONS

The viscosity of Glubran® 2 is only slightly higher than that of water; therefore, the product must be applied very carefully to prevent it from spreading to undesired areas. When necessary, apply gauze to protect the surrounding areas.

Always apply in minimal amounts, i.e. approximately one drop per square centimetre, and do not apply more than one drop to the same area. A second layer of Glubran® 2 may be applied over the first once it has polymerised.

Any excess product can be removed using a dry swab within 5–6 seconds of application.

The use of excessive amounts of Glubran® 2 in areas prone to infection, such as the vagina, has been associated with an increased risk of inflammation and/or infection. After solidification, any excess product may cause detachment of the adhesive layer and/or lead to the formation of small fragments that tend to detach from the tissues and must always be removed.

An excessive amount of product prolongs the setting time and may prevent adhesion.

Avoid contact with the eyes. In case of accidental contact, rinse immediately with water. If the product has polymerised, it will detach spontaneously after approximately 2–3 days.

Consult your doctor about possible side effects.

WARNINGS

-  Glubran® 2 must only be used by qualified doctors with experience in using the product. The manufacturer shall not be liable for damage caused by uses other than those indicated in this technical data sheet.
-  The product is for single use only.
-  Do not reuse the device after first opening. Reuse poses a substantial risk of infection to the patient due to the lack of sterility of the contents, failure of the adhesive properties and reduced effectiveness of the product due to contact with air.
-  The product is supplied ready for use.
-  The product must not be diluted or mixed with dyes or other substances, with the exception of triiodinated oil-based contrast media. Mixing with these substances proportionally affects polymerisation time; therefore, the mixture must be perfectly homogeneous and uniform.
-  Do not use Glubran® 2 with devices or accessories containing silicone or polycarbonate. Always ensure that the devices used are sterile and compatible with the product in order to prevent any induction of polymerisation or degradation of Glubran® 2.
-  Do not use the product if it is viscous and/or cloudy.
-  Any product residues must be disposed of.
-  The product cannot be re-sterilised.
-  Although rare, temporary local inflammatory reactions may occur after application.
-  When applied to the skin, although the exothermic reaction occurring during polymerisation does not exceed 45°C, Glubran® 2 may cause a slight burning sensation at the application site in particularly sensitive patients (e.g. children and the elderly).
-  Please refer to the instructions for use for the complete list of precautions and contraindications.

PRODUCT LIFESPAN

Following standard surgical procedures, the Glubran® 2 layer is eliminated over a period of time that varies depending on the type of tissue and the amount of Glubran® 2 applied. In embolisation procedures, Glubran® 2 persists for a longer period.

STERILITY

Glubran® 2 and its primary packaging are sterile and latex-free.

DISPOSAL

Once used, dispose of the device according to local procedures and guidelines.

CONTACTS

Any serious incident that occurs, related to the device, should be reported to the manufacturer GEM srl at info@gemitaly.it and to the competent authority of the member state where the user and/or patient is established. In this case, a detailed description of the problem occurred is required, as well as the variable data of the medical device involved (sales code, batch number, UDI code).

GEM S.r.l. Via dei Campi 2 - PO Box 427 - 55049 Viareggio (LU) Italy - Tel. +39 0584 389784/391388 - Fax +39 0584 397904
www.gemitaly.it - SH-GEM-info@hbfuller.com

SYMBOLS USED

Patient name

Please visit the website for further information

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